


New Foxy R1 and Foxy R2 Fraction Collectors







Launched at American Chemical Society, August 17-21 2008



Foxy R1: 68-3870-010 (replaces Foxy Jr)
Foxy R2: 68-2130-006 (replaces Foxy 200)


NEW Foxy racks!

- The Foxy R1 and Foxy R2 will ship with new racks, slightly lighter in color and RFID tags (for the R2).
- The Foxy R1 and Foxy R2 include choice of standard rack sizes. If no rack is specified, a 13mm test tube rack will be supplied.
- More expensive racks, i.e. prep racks, will have to be purchased as separate line items.

-  Automation
-  Flash Chromatography
-  Separations
-  Proteomics
-  Miniaturisation
-  E-Senses
-  Precision Fluid Delivery

...explore the possibilities

Tel: +44 (0) 1256 365492 Fax: +44 (0) 1256 365486
Email: enquiries@presearch.co.uk Web: www.presearch.co.uk



Presearch
A DIVISION OF UNITECH UK

Superior name-brand recognition with the following improvements:

- Sturdy cast aluminum construction
- Faster tube to tube transition time
- Practically zero void volume between the drop former and diverter valve
- Diverter valve, standard, for GLP requirements
- Simple touch screen user interface
- New rack for Foxy R2 holds (6) 96-well microplates
- Foxy R2 adds RFID technology for rack identification; fast and accurate setup
- Flow rates are 100 ml/min. A Foxy R2 version with flow rates up to 1000 ml/min and another Foxy R2 version with flow rates up to 5000 ml/min will be offered in the future

New Racks and Cables:

- All cables and racks currently provided for the Foxy 200 and Foxy Jr will be replaced with new cable and rack part numbers.
- Original Foxy racks can be used on the Foxy R1 and Foxy R2 without the R2-RFID advantage.
- A new rack for the Foxy R2 holds six 96-well microplates.

Literature for the Foxy R1 and Foxy R2 Fraction Collectors

- A datasheet is available.
- These fraction collectors and all of their accessories will be included on the September 1 Chromatography Price List.

Status of Foxy Jr and Foxy 200

Teledyne Isco will continue to build the Foxy Jr and Foxy 200 Fraction Collectors through October 2008. These fraction collectors will remain available while quantities last.

All stand alone Foxy Jr's and Foxy 200's will be removed from the September 1 Chromatography Price List.

Foxy 200 Systems

The following systems, including the Foxy 200 will be offered through December 2009. These systems will remain on the September 1 Chromatography Price List.

- Companion XL System
- Optix 10 System
- BioOptix 10 System



Automation



Flash Chromatography



Separations



Proteomics



Miniaturisation




E-Senses



Precision Fluid Delivery

...explore the possibilities

Tel: +44 (0) 1256 365492 Fax: +44 (0) 1256 365486
Email: enquiries@presearch.co.uk Web: www.presearch.co.uk



Presearch
A DIVISION OF UNITECH UK

Companion/Foxy Jr System

Companion Systems with the Foxy Jr will remain available as long as Foxy Jr's remain available. These Companion Systems will be removed from the September 1 Chromatography Price List.

Foxy Jr and Foxy 200 Systems with the Model 160 Gradient Formers

The new Foxy R1 and Foxy R2 fraction collectors do not support the Model 160 Gradient Formers. A few gradient formers are still available to sell.

- The ProTeam 210/320 Systems and the Programmable Density Gradient System including the Model 160 Gradient Former will be removed from the September 1 Chromatography Price List.

Status of Foxy Jr and Foxy 200 Accessories

Racks and cables

- The Foxy Jr and Foxy 200 racks will remain available until inventory is depleted. Once depleted, order entry will automatically substitute with new rack part number.
- Both Foxy Jr and Foxy 200 racks and cables will remain on the September 1 Chromatography Price List. All will be removed from the January 2009 price list.
- The Foxy Jr and Foxy 200 racks will be discontinued (obsolete) when inventory is depleted.
- The Foxy Jr and Foxy 200 cables will remain available for at least 5 years, not exceeding 7 years, based on availability of parts.

Diverter and Security valves








A diverter valve is included on ALL Foxy R1 and Foxy R2 Fraction Collectors.

- The Foxy Jr and Foxy 200 diverter and security valves will remain on the September 1 Chromatography Price List.
- These Foxy Jr and Foxy 200 valves will remain available for at least 5 years, not exceeding 7 years, based on availability of parts.

Model 160 Gradient Formers

The new Foxy R1 and R2 Fraction Collectors will NOT support our Model 160 Gradient Formers.

- The Model 160 Gradient Formers will be available only while quantities last.
- The Model 160 Gradient Formers will be removed from the September 1 Chromatography Price List.

-  Automation
-  Flash Chromatography
-  Separations
-  Proteomics
-  Miniaturisation
-  E-Senses
-  Precision Fluid Delivery

...explore the possibilities

Tel: +44 (0) 1256 365492 Fax: +44 (0) 1256 365486
Email: enquiries@presearch.co.uk Web: www.presearch.co.uk



Presearch
A DIVISION OF UNITECH UK